



Hemo Control n.x.t Generation microcuvettes

- ◆ Simple sampling from any angle
- ◆ Low 8 μ l sample volume
- ◆ Convenient package sizes
- ◆ Bubble-free technology

<p>Precise sample collection</p>	<ul style="list-style-type: none"> • Only 8 µl of capillary, venous or arterial blood • Pre-filled cuvette with dry reagent • Blood drawn in by capillary action • Minimised risk of air bubbles gives consistent results 	
<p>Individually packaged microcuvettes</p>	<ul style="list-style-type: none"> • Suitable for tropical conditions • Carton of 50. Ideal for mobile use and low volume sites • Each cuvette package contains its own desiccant • Low sensitivity to humidity • Operating temperature 15°C-40°C • Store at room temperature 15-30°C • Shelf life of 18 months 	
<p>Canister of 50 microcuvettes</p>	<ul style="list-style-type: none"> • Sealed canister with dessicant • Ideal for high volume sites • Operating temperature 15°C - 40°C • Store at room temperature 15-30°C • Shelf life of 24 months • Cuvettes stable for 3 months after opening when stored with the cap securely closed 	

Testing in three easy steps



1 Collect blood sample.



2 Put microcuvette into analyser.



3 Result appears in 25-60 seconds.

Distributed by

Manufacturer
 EKF-diagnostic GmbH
 Ebendorfer Chaussee 3
 39179 Barleben, Germany
 ☎ +49 (0)39203 5110
 sales@ekfdiagnostics.com

www.ekfdiagnostics.com

EKF | Diagnostics
 for life

Revision EN EU 4.0-01.2016